


Overdenture Attachments

Notch Implant GmbH collaborates with companies worldwide in order to offer its customers the highest quality available. By this approach, we supply our elastic caps for overdenture solutions from the world-renowned Italian company .

Retainers for Ball abutment



Color	Retention	Ref Code
Green	Elastic	AOD240140G
Yellow	Extra Soft	AOD240140Y
Pink	Soft	AOD240140P
Clear	Standart	AOD240140C
Black (Lab. Cap)	-	AOD240140B
Metal Housing	-	AOD240040M
Protective Disk	-	AOD240140D
Elastic Cap Pack	Multi	AOD240000P

Retainers for Locator abutment



Type	Color	Retention	Ref Code
With Pin	Yellow	Extra Soft	AOD260140Y
	Pink	Soft	AOD260140P
	Clear	Standart	AOD260140C
Without Pin	Violet	Hard	AOD260140V
	Yellow	Extra Soft	AOD250140Y
	Pink	Soft	AOD250140P
	Clear	Standart	AOD250140C
	Violet	Hard	AOD250140V
	Black (Lab. Cap)	-	AOD250140B
	Metal Housing	-	AOD250140M
Protective Disk	-	AOD250140D	
Elastic Cap Pack	Multi	AOD260000P	

Impression Coping

is used to make an impression of the positions of Locator abutments so that total denture is prepared in lab.

Lab Analog

is used to replicate positions of Locator abutments in stone cast.



System	Product	Ref Code
Locator	Impression Coping	PLI381010S
	Lab Analog	PLA381010S
Ball	Impression Coping	PBI381010S
	Lab Analog	PBA381010S